

## DEHYDRATION PROCEDURE LAYOUT

Generated By :ps

Date & Time : 27-02-2023 16:44:04

Sino	Protocol Name	Jar Name	Reagent/ Temperature	Incubation Time	No of Dips
1	DEHYDRATION PROCEDURE LAYOUT	J1	70%Alcohol	02:00	1
2	DEHYDRATION PROCEDURE LAYOUT	J2	95% Alcohol	02:00	1
3	DEHYDRATION PROCEDURE LAYOUT	J3	100 % Alcohol	02:00	1
4	DEHYDRATION PROCEDURE LAYOUT	J4	100 % Alcohol	02:00	1
5	DEHYDRATION PROCEDURE LAYOUT	J5	XYLENE-1	02:00	1
6	DEHYDRATION PROCEDURE LAYOUT	J6	XYLENE-2	02:00	1
7	DEHYDRATION PROCEDURE LAYOUT	J7	XYLENE-3	02:00	1